VOCO

Ionolux

Handmix
Instructions for use

Ionolux is a resin-reinforced light-curing radiopaque glass-polyalkenoate restorative. The material has outstanding aesthetics, good handling properties and high chemical adhesion.

Shades:
A1, A2, A3

Indications:

- Restorations class III and V, especially cervical fillings and root caries
- Fillings on deciduous teeth
- Small class I fillings
- Temporary fillings
- Core build-ups
- Linings

Method of Use:

Preparation:

In principle, cavity preparation should be carried out only minimally invasive, to save healthy tooth substance. It is not necessary to prepare undercuts. Class V cavities should have a minimum depth of 1 mm. Remove preparation residues with water spray, then dry with air so that the surface remains slightly glossy (do not overdry).

Lining

Areas close to the pulp should be protected with a calcium hydroxide lining.

Shade selection:

With the help of a VITA® shade guide on a moist tooth.

Mixing:

Shake the powder and the liquid bottle thoroughly for at least 3 s prior to dispensing. Dispense a level spoonful of powder onto a suitable mixing pad. Close the bottle immediately after dispensing to protect the powder against moisture. Dose 2 drops of liquid per level scoop of powder (powder/liquid ratio 3.2 : 1 g/g). Hold the bottle vertically when dispensing the drops. The dropper should be free from dried residues. Close the liquid bottle immediately after use. Work the first half of the powder into the liquid, then the second; for this, use spatulas of solid plastic or other abrasion-proof materials. Spread out the paste several times to achieve a homogeneous consistency.

Process the cement at room temperature (15°C – 23°C).

Application:

Place a matrix, if necessary. Apply the mixed **Ionolux** into the cavity with a suitable instrument and contour. Then light-cure with a suitable polymerisation lamp. Fillings of more than 2 mm thickness should be applied and light-cured in layers.

Light-curing:

The material can be light-cured with market-available polymerisation lamps. Their output should not be lower than 500 mW/cm². The curing time is 20 s per layer. Hold the tip of the light-curing lamp as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth will be smaller.

| | |
|-------------------------------|---------|
| Powder/liquid ratio [g/g] | 3.2:1.0 |
| Dosage [measuring spoon/drop] | 1:2 |
| Mixing time [s] | 30 |
| Working time [min] | >3:00 |
| Light-curing [s] | 20 |

Finishing:

After light-curing, remove excess material and finish and polish the filling as usual.

Information/Precautionary measures:

- In case of contact with oral mucosa or skin, remove immediately with a cotton soaked in alcohol. Flush with water.
- In case of contact with the eyes, rinse immediately with plenty of water and seek medical advice.
- Ionolux** contains polyacrylic acid, fluorosilicate glass, amines, BHT and methacrylates. Do not use **Ionolux** in case of known allergies against any of the ingredients.
- Phenolic substances, especially eugenol- and thymol-containing preparations disturb the polymerisation of the restorative. Do not use zinc-oxide eugenol cements or other eugenol-containing materials together with filling cements.

Storage:

Store at temperatures between 4°C – 23°C. Tightly close powder and liquid container immediately after use. Store dry and in an area protected from light.

Do not use **Ionolux** after the expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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