**VOCO**

# Merlon Plus QM

**GB****CE** 0482

In accordance with DIN EN ISO 9917-2

## Instructions for use

### Product description:

**Merlon Plus QM** is a radiopaque, self-curing, fluoride-releasing, resin-modified glass ionomer luting cement available in a paste-paste formula in the practical QuickMix syringe for quick and simple direct application. **Merlon Plus QM** offers optional removal of excess material following light-curing.

### Indications:

- Luting of
- metal-based inlays, onlays, crowns and bridges
  - high-strength (e.g., zirconium dioxide-based or lithium disilicates, etc.) all-ceramic inlays, onlays, crowns and bridges
  - metal, ceramic and glassfibre-reinforced posts
  - orthodontic bands and aids

### Method of use:

#### 1. Preparation:

Ensure that the surfaces to be cemented are clean. Thoroughly clean and dry enamel and dentine, and the inside surfaces of the restoration. Before luting, clean the stumps (e.g., pumice, water), rinse thoroughly and dry, ideally with a foam pellet or gentle blasts of air (take care not to overdry!). Prepare and treat the restoration as per the manufacturer's instructions.

Remove the protective cap from the QuickMix syringe. Due to technical reasons, there may be slight differences in the filling levels of the syringe chambers prior to first use. For this reason, when using a new syringe, dispense paste until the material exits both openings evenly. When reusing syringes, check the openings for any blockages and unblock as necessary. Attach a type 10 mixing tip in such a way that the notches on the mixing tip and the syringe are aligned. Turn the mixing tip 90° clockwise to lock the tip into place. To ensure even mixing of the two pastes, dispense and discard a little material before commencing the application.

#### Curing behaviour:

Working time:	2 min (23 °C, from start of mixing)
Removal of excess material:	a) in gel state (approx. 2 min after placement in mouth) or b) following light-curing (approx. 5 - 10 s per segment)
Setting time:	4.5 min (after placement in mouth)

#### 2. Placing the restoration:

Apply material to the preparation and/or the prepared restoration. Then place the restoration, applying gentle pressure so as to expel excess material. At 23 °C, the working time is 2 min from the start of mixing. Remove the used mixing tip and seal the syringe with the original twist cap until the next time it is required.

#### 3. Removal of excess material:

Excess cement is best removed in the gel phase, either following light-curing (5 - 10 s per segment) or following self-curing (approx. 2 min after placement), using a suitable instrument (e.g., scaler). Excess material in hard-to-reach areas should be removed first. Clean contact points and interproximal areas with dental floss. No material residue should remain in the sulcus, as this can result in irritation of the gingiva.

#### 4. Curing:

Curing is complete 4.5 min after placement. The restoration must be held in position during curing. Cured excess material can be removed with a finishing diamond or flexible polishing wheels. Finish the approximal areas with finishing/polishing strips.

### Information, precautionary measures:

**Merlon Plus QM** contains polyacrylic acid, peroxide, methacrylates and BHT. **Merlon Plus QM** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

In case of contact with the eyes, rinse immediately with plenty of clean water and seek medical advice.

Ensure sufficient residual dentine thickness to prevent pulp irritation. If necessary, protect the pulp or build up the core.

Do not use any phenolic substances (e.g., temporary cements containing eugenol), as these can impair the curing of **Merlon Plus QM**. Only use eugenol-free temporary cements for temporary luting.

Refit the original twist cap to prevent the **Merlon Plus QM** pastes drying out during storage.

The warranted product characteristics can only be achieved when using the original Sulzer Mixpac mixing tips as approved by **VOCO**. All the mixing tips for use with our respective products are available from **VOCO**.

### Storage:

Store **Merlon Plus QM** at temperatures of 4 °C to 23 °C with the original twist cap fitted. Do not use after expiry date.

Allow **Merlon Plus QM** to reach room temperature before use!

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB1727 E1 0118 99 © by VOCO
**VOCO GmbH**