PHILIPS



Sonic electric toothbrush

ProtectiveClean 4700

Built-in pressure sensor

2 cleaning modes

1 x BrushSync feature



HX6483/52



Cleaner teeth. Gently does it. Remove up to 10x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 10x more plaque than a manual toothbrush

Proven to improve oral health

Increased surface contact for eliminating plaque

Optimize your brushing

 \cdot Clean and gum care mode

Designed around your patient

- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly

Innovative technology

- Connects smart brush handle and smart brush heads
- A gentle nudge to go nice and easy
- · Patients gets the most from their brush heads
- · Philips Sonicare's advanced sonic technology

Removes up to 10x more plaque*



The C3 Premium Plaque Control brush head delivers our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact^{**} and removes up to 10x more plaque from hard-to-reach areas^{*}.

Too much pressure?



Patients may not notice if they are brushing too hard, but their toothbrush will. If patients start taking it out on their gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

Clean and gum care mode



Whether you want to focus on removing plaque or on getting your gums healthy and keeping them that way, this brush has a mode for you: Clean mode delivers superior cleaning, while Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums.

BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Smart replacement reminders



All brush heads wear out over time. But we've got that covered too: our smart brushrecognition technology tracks how long patients have been using their brush head for, and how hard they've been brushing. When it's time to replace it, a light on the handle lets them know. That way, they can be sure that their brush head is consistently doing a stellar job.

Easy-start program



Starting something new can take getting used to. That's why our easy-start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Helpful timers



Our QuadPacer lets patients know when they've spent just the right amount of time cleaning each part of their mouth, while our Smartimer tells them when they've brushed for the recommended two minutes.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Specifications

Modes

Clean: For exceptional everyday clean Gum Care: Gently massages gums

Items included

Handle: 1 ProtectiveClean Brush heads: 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W2c Optimal White compact Travel case: Hard travel case Charger: 1

Design and finishing

Color: White mint

Cleaning performance

Speed: Up to 62000 brush movement/min Performance: Removes up to 10x more plaque*

Health benefits: Helps reduce cavities Pressure feedback: Vibrates handle to alert user

Timer: Quadpacer and SmarTimer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

Ease of use

Handle compatibility: Easy click-on brush heads Battery indicator: Light shows battery status Handle: Slim ergonomic design Brushing time: Up to 2 weeks

Technical specifications

Battery: Rechargeable Battery type: Lithium ION Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

© 2021 Koninklijke Philips N.V. All Rights reserved.

their respective owners.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or

lssue date 2021-06-17 Version: 4.2.1

www.philips.com





* Removes up to 7x more plaque vs. a manual toothbrush * * compared to DiamondClean