

SMART|ONE ESSENTIAL FACE MASKS

Disposable 3-layer type IIR Face Masks with elastic band

CONFIGURATION OVERVIEW

- Designed for medical use
- Type IIR Face Mask
- Soft and comfortable
- Water repellent
- Fiberglass-free
- Latex-free
- PVC and Phthalates free
- Good performance of anti-micro-dust, anti-spray, anti-virus



Medical Device of Class I according to EU REG MDR 2017/745

CONFORM TO: EN ISO 10993-1:2021, EN ISO 10993-5:2009, EN ISO 10993-10:2010, ISO 15223-1:2021/Amd. 1:2025, EN ISO 14971:2019/A11:2021, EN ISO 20417:2021, EN 14683:2025

TECHNICAL SPECIFICATIONS:

- 3 layers = 2 polypropylene layers + 1 Melt blown layer (middle)
- Outer Layer: 25 g/m², PP, Coloured
- Middle Layer: 25 g/m², Purpose-made Meltblown, White
- Inner Layer: 25 g/m², PP nonwoven material, White
- 2 Holding Elastics: urethane elastic fiber
- Nose barret: Polyethylene with double iron wire bar

Dimensions of the mask (closed)	Dimensions of the mask (opened)	Length of nose barrette	Force at break for earloop
17.5 x 9.5 cm (± 0.5 cm)	17,5 x 15,5 cm (± 0.5 cm)	13 cm ±1 cm	1 Kg.

SHELF LIFE : 5 years

ARTICLE NUMBERS AND COLOUR VARIANTS:

FAMEA24-AKZ-BL-K02 = blue
 FAMEA24-AKZ-PI-K02 = pink
 FAMEA24-AKZ-WH-K02 = white



PACKAGING

Pieces/dispenser: 50
 Dispenser/carton: 40
 Carton size: 52 x 38 x 30 cm
 Carton weight: 8,10 Kg
 Cartons/layer : 4
 Layers/pallet : 6
 Cartons/pallet : 24
 Pallet size : 80x120x195 cm
 Pallet weight: 215 Kg



Akzenta International SA
 Via Giuseppe Motta, 24
 CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
 fax. +41 (0) 91 682 85 20

www.akzenta.com
 e-mail: business@akzenta.com

INSTRUCTIONS FOR USE:

1. Hold the mask and face the outer surface outward.
2. With the outer surface still facing away from you, place the ear loops comfortably around your ears.
3. Adjust the noseband to fit securely over the bridge of your nose, ensuring it conforms to the shape of your face.
4. While holding the noseband in place, pull the bottom edge of the mask under your chin, making sure it fully covers your nose, mouth, and chin.



Rev. 3 03-07-2026



Akzenta International SA
Via Giuseppe Motta, 24
CH-6830 Chiasso

tel. +41 (0) 91 921 14 92
fax. +41 (0) 91 682 85 20

www.akzenta.com
e-mail: business@akzenta.com